

Product Datasheet

InVivoMAb Anti-SARS-CoV-2 NP/Nucleoprotein Antibody (orb1290135)

Catalog Number	orb1290135
Category	Antibodies
Description	Human monoclonal antibody to SARS-CoV Nucleocapsid
Clonality	Recombinant
Isotype	Human IgG
Conjugation	Unconjugated
Form/Appearance	Liquid
Buffer/Preservatives	PBS,pH7.5
Immunogen	Recombinant SARS-CoV Nucleocapsid protein
Dilution range	WB:1:1000-2000; ELISA:1:5000-1:10000
Specificity	Recognizes SARS-CoV Nucleocapsid protein
Antibody Type	Recombinant Antibody
Clone Number	BUE25X221XET5221
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)