

Product Datasheet

Balstilimab (orb1290180)

Catalog Number	orb1290180
Category	Antibodies
Description	Balstilimab biosimilar antibody
Target	PDCD1/PD1/CD279
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG4-kappa
Conjugation	Unconjugated
Reactivity	Human
Buffer/Preservatives	PBS buffer pH 7.5
Purity	> 95%
Tested applications	ELISA, FC
Antibody Type	Biosimilar Antibody
Source	CHO cells
CAS Number	2230167-06-1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558