

Product Datasheet

MPXV A29L Rabbit Polyclonal Antibody (orb1291565)

Catalog Number	orb1291565
Category	Antibodies
Description	MPXV A29L Rabbit Polyclonal Antibody
Target	M1R
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Virus
Predicted Reactivity	Virus
Form/Appearance	Liquid
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	Recombinant MPXV/A29L protein (1-110/110aa)
UniProt ID	Q77HM6
MW	12.5 kDa

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558


Tested applications	ELISA
Dilution range	ELISA=1:5000-10000
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Rabbit Anti-MPXV A29L Antibody Bind with MPXV A29L protein, His Tag



Measured by its binding ability in a indirect ELISA. Immobilized MPXV A29L protein, His Tag at 2 µg/mL (100 µl/well) can bind Rabbit Anti-MPXV A29L Antibody, the EC50 is 150.9 ng/ml.

Rabbit Anti-MPXV A29L Antibody Bind with MPXV A29L protein, His Tag



Measured by its binding ability in a indirect ELISA. Immobilized MPXV A29L protein, His Tag at 2 µg/mL (100 µl/well) can bind Rabbit Anti-MPXV A29L Antibody, the EC50 is 42.93 ng/ml.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558