

Product Datasheet

L-asarinin (orb1303777)

Description	L-asarinin
Storage	-20°C
Note	For research use only
Formula	$C_{20}H_{18}O_6$
Purity	98%
MW	354.35
CAS Number	133-05-1
SMILES	<chem>C1Oc2ccc(cc2O1)[C@H]1OC[C@H]2[C@@H]1CO[C@H]2c1ccc2OCOc2c1</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)