

Product Datasheet

S Protein Antibody (orb1312103)

Catalog Number	orb1312103
Category	Antibodies
Description	Mouse monoclonal antibody to S Protein
Target	S Protein
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2b
Conjugation	Unconjugated
Reactivity	Virus
Immunogen	Recombinant SARS-CoV-2 Spike-RBD protein produced in E.coli.
MW	34kDa
Tested applications	ELISA
Dilution range	ELISA 1:5000-10000
Antibody Type	Primary Antibody
Clone Number	OTI12F9
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	YP_009724390

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)