

Product Datasheet

Purified HIV1 p24 Mouse Monoclonal Antibody [Clone ID: OTI1A1] (orb1325671)

Catalog Number	orb1325671
Category	Antibodies
Description	Mouse monoclonal antibody to Purified HIV1 p24
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2a
Conjugation	Unconjugated
Reactivity	Virus
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 1-211 of HIV gag p24
Tested applications	WB
Dilution range	WB: 1:500-1:2000
Antibody Type	Primary Antibody
Clone Number	OTI1A1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)