

## Product Datasheet

### VWA5A Antibody (orb1337404)

|                             |   |
|-----------------------------|---|
| <b>Catalog Number</b>       | orb1337404  |
| <b>Category</b>             | Antibodies  |
| <b>Description</b>          | VWA5A Antibody  |
| <b>Target</b>               | VWA5A   |
| <b>Clonality</b>            | Monoclonal  |
| <b>Species/Host</b>         | Mouse   |
| <b>Isotype</b>              | IgG1  |
| <b>Conjugation</b>          | Unconjugated  |
| <b>Reactivity</b>           | Human, Monkey, Mouse  |
| <b>Buffer/Preservatives</b> | Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)   |
| <b>Purification</b>         | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)                               |
| <b>Immunogen</b>            | Full length human recombinant protein of human VWA5A(NP_055437) produced in HEK293T cell.   |
| <b>Tested applications</b>  | FC, IF, WB  |
| <b>Dilution range</b>       | WB 1:500~1000, IF 1:100, FLOW 1:100   |
| <b>Clone Number</b>         | BQVK5C5   |
| <b>Storage</b>              | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**


68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Note** For research use only

**Expiration Date** 12 months from date of receipt.



Western blot analysis of HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY VWA5A (Right lane) cDNA for 48 hrs and lysed using Breast cancer suppressor candidate 1 (VWA5A) antibody

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)