

Product Datasheet

LTA Human Over-expression Lysate (orb1339694)

| | |
|--------------------------|--|
| Catalog Number | orb1339694 |
| Category | Proteins |
| Description | Transient overexpression lysate of lymphotoxin alpha (TNF superfamily, member 1) (LTA), transcript variant 1 |
| Target | LTA |
| Tag | C-Myc/DDK |
| UniProt ID | P01374 |
| MW | 22.3 kDa |
| Expression System | HEK293T |
| Source | Human |
| Storage | Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. |
| Dry Ice Shipping | Please note: This product requires shipment on dry ice. A dry ice surcharge will apply. |
| Note | For research use only |
| NCBI | NM_001159740, NP_001153212 |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)