

## Product Datasheet

### p73 (TP73) Human Over-expression Lysate (orb1341217)

|                          |                                                                                                                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb1341217                                                                                                                                                                                                                                               |
| <b>Category</b>          | Proteins                                                                                                                                                                                                                                                 |
| <b>Description</b>       | Transient overexpression lysate of tumor protein p73 (TP73), transcript variant 4                                                                                                                                                                        |
| <b>Target</b>            | p73                                                                                                                                                                                                                                                      |
| <b>Tag</b>               | C-Myc/DDK                                                                                                                                                                                                                                                |
| <b>UniProt ID</b>        | <b>O15350</b>                                                                                                                                                                                                                                            |
| <b>MW</b>                | 47.5 kDa                                                                                                                                                                                                                                                 |
| <b>Expression System</b> | HEK293T                                                                                                                                                                                                                                                  |
| <b>Source</b>            | Human                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. |
| <b>Dry Ice Shipping</b>  | <b>Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.</b>                                                                                                                                                           |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                    |
| <b>NCBI</b>              | <b>NP_001119714, NM_001126242</b>                                                                                                                                                                                                                        |

#### Biorbyt Ltd.


7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Western blot analysis of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using p73 (TP73) Human Over-expression Lysate

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)