


## Product Datasheet

# Integrin alpha 3 (ITGA3) Human Over-expression Lysate (orb1346557)

|                        |                                                                                                                                                                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Transient overexpression lysate of integrin, alpha 3 (antigen CD49C, alpha 3 subunit of VLA-3 receptor) (ITGA3), transcript variant b                                                                                                                    |
| <b>Target</b>          | Integrin alpha 3                                                                                                                                                                                                                                         |
| <b>Storage</b>         | Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. |
| <b>Tag</b>             | C-Myc/DDK                                                                                                                                                                                                                                                |
| <b>Note</b>            | For research use only                                                                                                                                                                                                                                    |
| <b>Source</b>          | Human                                                                                                                                                                                                                                                    |
| <b>MW</b>              | 118.79 kDa                                                                                                                                                                                                                                               |
| <b>Uniprot ID</b>      | <b>P26006</b>                                                                                                                                                                                                                                            |
| <b>NCBI</b>            | <b>NM_005501, NP_005492</b>                                                                                                                                                                                                                              |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                                                                                                                                           |



Western blot analysis of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using Integrin alpha 3 (ITGA3) Human Over-expression Lysate

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)