

Product Datasheet

AIFM3 Human Over-expression Lysate (orb1352918)

Catalog Number	orb1352918
Category	Proteins
Description	Transient overexpression lysate of apoptosis-inducing factor, mitochondrion-associated, 3 (AIFM3), nuclear gene encoding mitochondrial protein, transcript variant 1
Target	AIFM3
Tag	C-Myc/DDK
UniProt ID	Q96NN9
MW	66.6 kDa
Expression System	HEK293T
Source	Human
Storage	Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
Dry Ice Shipping	Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.
Note	For research use only
NCBI	NM_144704, NP_653305
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)