

## Product Datasheet

### DEFB106A Human Over-expression Lysate (orb1353271)

|                   |                                                                                                                                                                                                                                                          |
|-------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number    | orb1353271                                                                                                                                                                                                                                               |
| Category          | Proteins                                                                                                                                                                                                                                                 |
| Description       | Transient overexpression lysate of defensin, beta 106A (DEFB106A)                                                                                                                                                                                        |
| Target            | DEFB106A                                                                                                                                                                                                                                                 |
| Tag               | C-Myc/DDK                                                                                                                                                                                                                                                |
| UniProt ID        | <b>Q8N104</b>                                                                                                                                                                                                                                            |
| MW                | 7.4 kDa                                                                                                                                                                                                                                                  |
| Expression System | HEK293T                                                                                                                                                                                                                                                  |
| Source            | Human                                                                                                                                                                                                                                                    |
| Storage           | Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. |
| Dry Ice Shipping  | <b>Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.</b>                                                                                                                                                           |
| Note              | For research use only                                                                                                                                                                                                                                    |
| NCBI              | <b>NM_152251, NP_689464</b>                                                                                                                                                                                                                              |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)