

## Product Datasheet

# TAFA4 Human Over-expression Lysate (orb1355800)

|                          |  |
|--------------------------|--|
| <b>Catalog Number</b>    | orb1355800   |
| <b>Category</b>          | Proteins   |
| <b>Description</b>       | Transient overexpression lysate of family with sequence similarity 19 (chemokine (C-C motif)-like), member A4 (FAM19A4), transcript variant 1  |
| <b>Target</b>            | TAFA4  |
| <b>Tag</b>               | C-Myc/DDK  |
| <b>UniProt ID</b>        | <b>Q96LR4</b>  |
| <b>MW</b>                | 15.7 kDa   |
| <b>Expression System</b> | HEK293T  |
| <b>Source</b>            | Human  |
| <b>Storage</b>           | Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C. |
| <b>Dry Ice Shipping</b>  | <b>Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.</b>   |
| <b>Note</b>              | For research use only  |
| <b>NCBI</b>              | <b>NP_872328, NM_182522</b>  |
| <b>Expiration Date</b>   | 6 months from date of receipt.   |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)