

## Product Datasheet

# SLC34A2 Human Over-expression Lysate (orb1357328)

<b>Catalog Number</b>	orb1357328
<b>Category</b>	Proteins
<b>Description</b>	Transient overexpression lysate of solute carrier family 34 (sodium phosphate), member 2 (SLC34A2)
<b>Target</b>	SLC34A2
<b>Tag</b>	C-Myc/DDK
<b>UniProt ID</b>	<b>O95436</b>
<b>MW</b>	75.5 kDa
<b>Expression System</b>	HEK293T
<b>Source</b>	Human
<b>Storage</b>	Ship with dry ice. Upon receiving, store the sample at -80°C. Avoid repeated freeze-thaw cycles. Lysate samples can be diluted with 2xSDS Sample Buffer provided. Lysate samples are stable for 12 months from the date of receipt when stored at -80°C.
<b>Dry Ice Shipping</b>	<b>Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.</b>
<b>Note</b>	For research use only
<b>NCBI</b>	<b>NP_006415, NM_006424</b>
<b>Expiration Date</b>	6 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)