

## Product Datasheet

# PSMA (FOLH1) Human Over-expression Lysate (orb1359185)

|                          |                                                                                                                                                                                                                                                                                |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb1359185                                                                                                                                                                                                                                                                     |
| <b>Category</b>          | Proteins                                                                                                                                                                                                                                                                       |
| <b>Description</b>       | FOLH1 HEK293T cell transient overexpression lysate (as WB positive control)                                                                                                                                                                                                    |
| <b>Target</b>            | PSMA                                                                                                                                                                                                                                                                           |
| <b>Tag</b>               | C-Myc/DDK                                                                                                                                                                                                                                                                      |
| <b>UniProt ID</b>        | <b>Q04609</b>                                                                                                                                                                                                                                                                  |
| <b>MW</b>                | 50.5 kDa                                                                                                                                                                                                                                                                       |
| <b>Expression System</b> | HEK293T                                                                                                                                                                                                                                                                        |
| <b>Source</b>            | Human                                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | Ship at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C. |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                          |
| <b>NCBI</b>              | <b>NM_001193473, NP_001180402</b>                                                                                                                                                                                                                                              |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                                                                                 |

### Biorbyt Ltd.


7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Western blot validation of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using PSMA (FOLH1) Human Over-expression Lysate

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)