

## Product Datasheet

### LOC136242 (PRSS37) Human Over-expression Lysate (orb1362824)


|                        |                                                                                                                                                                                                                                                                                |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>  | orb1362824                                                                                                                                                                                                                                                                     |
| <b>Description</b>     | PRSS37 HEK293T cell transient overexpression lysate (as WB positive control)                                                                                                                                                                                                   |
| <b>Target</b>          | LOC136242                                                                                                                                                                                                                                                                      |
| <b>Storage</b>         | Ship at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C. |
| <b>Tag</b>             | C-Myc/DDK                                                                                                                                                                                                                                                                      |
| <b>Note</b>            | For research use only                                                                                                                                                                                                                                                          |
| <b>Source</b>          | Human                                                                                                                                                                                                                                                                          |
| <b>MW</b>              | 26.3 kDa                                                                                                                                                                                                                                                                       |
| <b>Uniprot ID</b>      | <b>A4D1T9</b>                                                                                                                                                                                                                                                                  |
| <b>NCBI</b>            | <b>NM_001008270, NP_001008271</b>                                                                                                                                                                                                                                              |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                                                                                                                                                                 |

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Western blot validation of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using LOC136242 (PRSS37) Human Over-expression Lysate

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)