

## Product Datasheet

### IL32 Human Over-expression Lysate (orb1363180)

|                          |                                                                                                                                                                                                                                                                                |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb1363180                                                                                                                                                                                                                                                                     |
| <b>Category</b>          | Proteins                                                                                                                                                                                                                                                                       |
| <b>Description</b>       | IL32 HEK293T cell transient overexpression lysate (as WB positive control)                                                                                                                                                                                                     |
| <b>Target</b>            | IL32                                                                                                                                                                                                                                                                           |
| <b>Tag</b>               | C-Myc/DDK                                                                                                                                                                                                                                                                      |
| <b>UniProt ID</b>        | <b>P24001</b>                                                                                                                                                                                                                                                                  |
| <b>MW</b>                | 21.7 kDa                                                                                                                                                                                                                                                                       |
| <b>Expression System</b> | HEK293T                                                                                                                                                                                                                                                                        |
| <b>Source</b>            | Human                                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | Ship at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C. |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                                                          |
| <b>NCBI</b>              | <b>NP_001012736, NM_001012718</b>                                                                                                                                                                                                                                              |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                                                                                 |

#### Biorbyt Ltd.


7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Western blot validation of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using IL32 Human Over-expression Lysate

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)