


## Product Datasheet

# Uroguanylin (GUCA2B) Human Over-expression Lysate (orb1367675)

|                        |                                                                                                                                                                                                                                                                                |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | GUCA2B HEK293T cell transient overexpression lysate (as WB positive control)                                                                                                                                                                                                   |
| <b>Target</b>          | Uroguanylin                                                                                                                                                                                                                                                                    |
| <b>Storage</b>         | Ship at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C. |
| <b>Tag</b>             | C-Myc/DDK                                                                                                                                                                                                                                                                      |
| <b>Note</b>            | For research use only                                                                                                                                                                                                                                                          |
| <b>Source</b>          | Human                                                                                                                                                                                                                                                                          |
| <b>MW</b>              | 12.1 kDa                                                                                                                                                                                                                                                                       |
| <b>Uniprot ID</b>      | <a href="#">Q16661</a>                                                                                                                                                                                                                                                         |
| <b>NCBI</b>            | <a href="#">NP_009033</a> , <a href="#">NM_007102</a>                                                                                                                                                                                                                          |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                                                                                                                                                                                                                 |



Western blot validation of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using Uroguanylin (GUCA2B) Human Over-expression Lysate

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)