

## Product Datasheet

# MARCHF2 Human Over-expression Lysate (orb1369283)

|                   |                                                                                                                                                                                                                                                                                |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog Number    | orb1369283                                                                                                                                                                                                                                                                     |
| Category          | Proteins                                                                                                                                                                                                                                                                       |
| Description       | 42065 HEK293T cell transient overexpression lysate (as WB positive control)                                                                                                                                                                                                    |
| Target            | MARCHF2                                                                                                                                                                                                                                                                        |
| Tag               | C-Myc/DDK                                                                                                                                                                                                                                                                      |
| UniProt ID        | <b>Q9P0N8</b>                                                                                                                                                                                                                                                                  |
| MW                | 27 kDa                                                                                                                                                                                                                                                                         |
| Expression System | HEK293T                                                                                                                                                                                                                                                                        |
| Source            | Human                                                                                                                                                                                                                                                                          |
| Storage           | Ship at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C. |
| Note              | For research use only                                                                                                                                                                                                                                                          |
| NCBI              | <b>NP_057580, NM_016496</b>                                                                                                                                                                                                                                                    |
| Expiration Date   | 6 months from date of receipt.                                                                                                                                                                                                                                                 |

### Biorbyt Ltd.


7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Western blot validation of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using MARCHF2 Human Over-expression Lysate

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)