

Product Datasheet

Sigma1 receptor (SIGMAR1) Human Over-expression Lysate (orb1373511)

Catalog Number	orb1373511
Category	Proteins
Description	SIGMAR1 HEK293T cell transient overexpression lysate (as WB positive control)
Target	Sigma1 receptor
Tag	C-Myc/DDK
UniProt ID	Q99720
MW	21.3 kDa
Expression System	HEK293T
Source	Human
Storage	Ship at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
Note	For research use only
NCBI	NM_147157, NP_671513
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)