

Product Datasheet

Kv beta 2 (KCNAB2) Human Over-expression Lysate (orb1374172)

Catalog Number	orb1374172
Category	Proteins
Description	KCNAB2 HEK293T cell transient overexpression lysate (as WB positive control)
Target	Kv beta 2
Tag	C-Myc/DDK
UniProt ID	Q13303
MW	39.1 kDa
Expression System	HEK293T
Source	Human
Storage	Ship at ambient temperature. Upon receiving, store the sample at -20°C. Lysate samples can be reconstituted with SDS Sample Buffer. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples are stable for 12 months from date of receipt when stored at -20°C.
Note	For research use only
NCBI	NP_742128, NM_172130
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Western blot validation of Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with transfection reagent MegaTran 2.0) using Kv beta 2 (KCNAB2) Human Over-expression Lysate

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)