

## Product Datasheet

### Human TNFRSF17/BCMA ELISA Kit (orb138060)

|                            |                                                                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb138060                                                                                                                                                                         |
| <b>Category</b>            | Assays and Kits                                                                                                                                                                   |
| <b>Description</b>         | Human TNFRSF17/BCMA ELISA Kit                                                                                                                                                     |
| <b>Target</b>              | Tumor necrosis factor receptor superfamily member 17                                                                                                                              |
| <b>Reactivity</b>          | Human                                                                                                                                                                             |
| <b>Immunogen</b>           | Expression system for standard: NS0; Immunogen sequence: M1-A54                                                                                                                   |
| <b>UniProt ID</b>          | <b>Q02223</b>                                                                                                                                                                     |
| <b>Tested applications</b> | ELISA                                                                                                                                                                             |
| <b>Specificity</b>         | Natural and recombinant human TNFRSF17                                                                                                                                            |
| <b>Cross Reactivity</b>    | There is no detectable cross-reactivity with other relevant proteins.                                                                                                             |
| <b>Range</b>               | 31.2 pg/ml - 2,000 pg/ml                                                                                                                                                          |
| <b>Sensitivity</b>         | 10 pg/ml                                                                                                                                                                          |
| <b>Sample Types</b>        | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).                                                                                                        |
| <b>Storage</b>             | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| <b>Note</b>                | For research use only                                                                                                                                                             |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                    |

#### Biorbyt Ltd.


7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Human TNFRSF17 ELISA Kit**

Human TNFRSF17/BCMA PicoKine ELISA Kit standard curve

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)