

## Product Datasheet

### EIF4ENIF1 Antibody (N-term) (orb1436565)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb1436565                                                                                                                                                                   |
| <b>Category</b>            | Antibodies                                                                                                                                                                   |
| <b>Description</b>         | Rabbit polyclonal antibody to EIF4ENIF1.                                                                                                                                     |
| <b>Target</b>              | EIF4ENIF1                                                                                                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                                                                        |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                            |
| <b>UniProt ID</b>          | <b>Q9NRA8</b>                                                                                                                                                                |
| <b>MW</b>                  | 108201                                                                                                                                                                       |
| <b>Tested applications</b> | WB                                                                                                                                                                           |
| <b>Dilution range</b>      | WB - 1:1000                                                                                                                                                                  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles                                       |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

#### Biorbyt LLC


68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

**Note** For research use only

**NCBI** **NP\_001157974.1, NP\_001157973.1, NP\_062817.2**

**Expiration Date** 12 months from date of receipt.



EIF4ENIF1 Antibody (N-term) western blot analysis in WiDr cell line lysates (35 ug/lane). This demonstrates the EIF4ENIF1 antibody detected the EIF4ENIF1 protein (arrow).

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)