

## Product Datasheet

### Kallikrein 9 Antibody Blocking peptide (orb1448049)

|                          |                                                                                                                                                                                                                                                     |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb1448049                                                                                                                                                                                                                                          |
| <b>Category</b>          | Proteins                                                                                                                                                                                                                                            |
| <b>Description</b>       | Kallikrein 9 Antibody Blocking peptide                                                                                                                                                                                                              |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.                                                                                                                |
| <b>Protein Sequence</b>  | The synthetic peptide sequence used to generate the antibody orb1930082 was selected from the N-term region of human KLK9. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| <b>UniProt ID</b>        | <b>Q9UKQ9</b>                                                                                                                                                                                                                                       |
| <b>Application notes</b> | Synthetic peptide                                                                                                                                                                                                                                   |
| <b>Note</b>              | For research use only                                                                                                                                                                                                                               |

---

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)