

Product Datasheet

BRAF (Phospho-S445) Antibody Blocking peptide (orb1450045)

Catalog Number	orb1450045
Category	Proteins
Description	BRAF (Phospho-S445) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1931167 was selected from the region of human Phospho-BRAF-pS445. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
UniProt ID	P15056
Application notes	Synthetic peptide
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)