

Product Datasheet

Bi-Syk (Phospho-Y525/526) Antibody Blocking peptide (orb1450174)

Catalog Number	orb1450174
Category	Proteins
Description	Bi-Syk (Phospho-Y525/526) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Protein Sequence	The synthetic peptide sequence was selected from the region of human Phospho-Syk-Y525/526. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
UniProt ID	P43405
Application notes	Synthetic peptide
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)