

Product Datasheet

Bi-FGFR(Phospho-Y653/654) Antibody Blocking peptide (orb1450290)

Catalog Number	orb1450290
Category	Proteins
Description	Bi-FGFR(Phospho-Y653/654) Antibody Blocking peptide
Target	FGFR1
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Protein Sequence	The synthetic peptide sequence was selected from the 649~658 region of human Phospho-FGFR-Y653/Y654. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
UniProt ID	P11362
Application notes	Synthetic peptide
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)