

Product Datasheet

Cyclin D1 (S90) Antibody Blocking peptide (orb1450779)

Catalog Number	orb1450779
Category	Proteins
Description	Cyclin D1 (S90) Antibody Blocking peptide
Target	CCND1 {ECO:0000303 PubMed:8204893, ECO:0000312 HGNC:HGNC:1582}
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1931553 was selected from the S90 region of human Cyclin D1. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
UniProt ID	P24385
Application notes	Synthetic peptide
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)