

Product Datasheet

LC3 (APG8C) (T48) Antibody Blocking peptide (orb1453163)

Catalog Number	orb1453163
Category	Proteins
Description	LC3 (APG8C) (T48) Antibody Blocking peptide
Target	MAP1LC3C
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1933536 was selected from the MAP1LC3C region of human APG8c (MAP1LC3C). A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
UniProt ID	Q9BXW4
Application notes	Synthetic peptide
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)