

Product Datasheet

LC3 (APG8A) (P45) Antibody Blocking peptide (orb1453209)

| | |
|--------------------------|--|
| Catalog Number | orb1453209 |
| Category | Proteins |
| Description | LC3 (APG8A) (P45) Antibody Blocking peptide |
| Target | MAP1LC3A |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Protein Sequence | The synthetic peptide sequence used to generate the antibody orb1933560 was selected from the region of human LC3 (APG8a) beginning with residue P45. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| UniProt ID | Q9H492 |
| Application notes | Synthetic peptide |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)