

## **Product Datasheet**

## HSP70/HSC70 Antibody: HRP (orb146754)

Catalog Number orb146754

**Category** Antibodies

**Description** Mouse monoclonal to Hsp70 (HRP). Hsp70 genes encode abundant heat-

inducible 70-kDa hsps (hsp70s). In most eukaryotes hsp70 genes exist as part of a multigene family. They are found in most cellular compartments of eukaryotes including nuclei, mitochondria, chloroplasts, the endoplasmic reticulum and the cytosol, as well as in bacteria. The genes show a high degree of conservation, having at least 50% identity. The N-terminal two thirds of hsp70s are more conserved than the C-terminal third. Hsp70 binds ATP with high affinity and possesses a weak ATPase activity which can be stimulated by binding to unfolded proteins and synthetic peptides. When hsc70 (constitutively expressed) present in mammalian cells was truncated, ATP binding activity was found to reside in an N-terminal fragment of 44 kDa which lacked peptide binding capacity. Polypeptide binding ability therefore resided within the C-terminal half. The structure of this ATP binding domain displays multiple features of nucleotide binding proteins. All hsp70s, regardless of location, bind proteins, particularly unfolded ones. The molecular chaperones of the hsp70 family recognize and bind to nascent polypeptide chains as well as partially folded intermediates of proteins preventing their aggregation and misfolding. The binding of ATP triggers a critical conformational change leading to the release of the bound substrate protein. The universal ability of hsp70s to undergo cycles of binding to and release from hydrophobic stretches of partially unfolded proteins determines their role in a great variety of vital intracellular functions such as protein synthesis, protein folding and oligomerization and protein transport.

Target HSP70/HSC70

**Clonality** Monoclonal

Species/Host Mouse

**Isotype** lgG1

**Conjugation** HRP



## **Biorbyt.com**

**Reactivity** Bovine, Canine, Drosophila, Fish, Frog, Gallus, Guinea pig, Hamster, Human,

Monkey, Mouse, Other, Plant, Porcine, Rabbit, Rat, Sheep

**Concentration** 1 mg/ml

**Buffer/Preservatives** 73.64mM Carbonate, 54.55mM Ethanolamine, 45.45mM Cyanoborohydride,

18.18mM Sodium Hydroxide and 0.23mM Citrate in dH2O

**Purification** Protein G Purified

Immunogen Recombinant HSP70/HSC70

UniProt ID PODMV9, PODMV8

MW 73kDa

**Tested applications** ELISA, IHC, WB

**Dilution range** WB (1:1000), IHC (1:100), ICC/IF (1:50)

**Application notes** 1 μg/ml of SMC-104 was sufficient for detection of HSP70/HSC70 in 20 μg of heat

shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-

mouse IgG:HRP as the secondary antibody.

Clone Number N27F3-4

**Storage** Conjugated antibodies should be stored according to the product label

**Note** For research use only

Entrez 3303


NCBI NP\_005336.3

**Expiration Date** 12 months from date of receipt.

Phone: <u>+1 (415) 906-5211</u> | Fax: <u>+1 (415) 651-8558</u>







Immunohistochemistry analysis using Mouse Anti-Hsp70 Monoclonal Antibody, Clone N27F3-4. Tissue: backskin. Species: Mouse. Fixation: Bouin's Fixative and paraffin-embedded. Primary Antibody: Mouse Anti-Hsp70 Monoclonal Antibody at 1:100 for 1 hour at RT. Secondary Antibody: FITC Goat Anti-Mouse (green) at 1:50 for 1 hour at RT. Localization: Epidermis.



Western Blot analysis of Human Cell lysates showing detection of Hsp70 protein using Mouse Anti-Hsp70 Monoclonal Antibody, Clone N27F3-4. Load: 15  $\mu$ g. Block: 1.5% BSA for 30 minutes at RT. Primary Antibody: Mouse Anti-Hsp70 Monoclonal Antibody at 1:1000 for 2 hours at RT. Secondary Antibody: Sheep Anti-Mouse IgG: HRP for 1 hour at RT.

Phone: <u>+1 (415) 906-5211</u> | Fax: <u>+1 (415) 651-8558</u>