


Product Datasheet

Myelin Protein Zero Rabbit mAb (orb1472217)

Catalog Number	orb1472217
Description	Myelin Protein Zero Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Uniprot ID	P25189
Expiration Date	12 months from date of receipt.



Western blot detection of Myelin Protein Zero in C6 cell lysates using Myelin Protein Zero Rabbit mAb (1:1000 diluted).
Predicted band size: 28 kDa. Observed band size: 28 kDa.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)