

Product Datasheet

Goat Mouse IgM - Texas Red Antibody (orb1473571)

Catalog Number	orb1473571
Category	Antibodies
Description	The Goat Mouse IgM - Texas Red Antibody is suitable for ELISA, FC, IF, IHC. It is a Polyclonal, Texas Red conjugated antibody which raised against Mouse IgM. Purification: The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Clonality	Polyclonal
Species/Host	Goat
Conjugation	Texas Red
Reactivity	Mouse
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Immunogen	Mouse IgM
Tested applications	ELISA, FC, IF, IHC
Dilution range	IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000)
Antibody Type	Secondary Antibody
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)