

Product Datasheet

JNK1 Rabbit Monoclonal Antibody (orb1474519)

Catalog Number	orb1474519
Category	Antibodies
Description	The JNK1 Antibody is suitable for IF, IP, WB. It is a Monoclonal, Unconjugated antibody which raised against Recombinant protein of human JNK1 .Purification: The antibody was purified by immunogen affinity chromatography.
Target	MAPK8
Clonality	Monoclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Mouse, Rat
Form/Appearance	Liquid
Buffer/Preservatives	50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin.
Purification	The antibody was purified by immunogen affinity chromatography.
Immunogen	Recombinant protein of human JNK1
UniProt ID	Q91Y86, P49185
Tested applications	IF, IP, WB
Dilution range	WB (1/500 - 1/1000), IF/IC (1/50 - 1/100), IP (1/10 - 1/50)
Specificity	Recognizes endogenous levels of JNK1 protein.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Entrez	26419, 116554
Expiration Date	12 months from date of receipt.



Western blot analysis of JNK1 expression in C6 (A), NIH3T3 (B) whole cell lysates. (Predicted band size: 48 kD; Observed band size: 46, 54 kD)



Immunofluorescent analysis of JNK1 staining in HeLa cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a AF488-conjugated secondary antibody (green) in PBS at room temperature in the dark.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)