

Product Datasheet

Importin 9 Rabbit Monoclonal Antibody (orb1474624)

| | |
|-----------------------------|---|
| Catalog Number | orb1474624 |
| Category | Antibodies |
| Description | The Importin 9 Antibody is suitable for IF, IP, WB. It is a Monoclonal, Unconjugated antibody which raised against A synthetic peptide of human Importin 9 .Purification: The antibody was purified by immunogen affinity chromatography. |
| Target | IPO9 |
| Clonality | Monoclonal |
| Species/Host | Rabbit |
| Conjugation | Unconjugated |
| Reactivity | Human |
| Form/Appearance | Liquid |
| Buffer/Preservatives | 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin. |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Immunogen | A synthetic peptide of human Importin 9 |
| UniProt ID | Q96P70 |
| Tested applications | IF, IP, WB |
| Dilution range | WB (1/500 - 1/1000), IF/IC (1/50 - 1/100), IP (1/10 - 1/50) |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Specificity | Recognizes endogenous levels of Importin 9 protein. |
| Antibody Type | Primary Antibody |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Entrez | 55705 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Importin 9 expression in Jurkat (A), HeLa (B) whole cell lysates. (Predicted band size: 116 kD; Observed band size: 116 kD)



Immunofluorescent analysis of Importin 9 staining in A549 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a AF488-conjugated secondary antibody (green) in PBS at room temperature in the dark.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)