

Product Datasheet

Rabbit Human IgG (H+L) Secondary Antibody, FITC Conjugated (orb1474920)

Catalog Number	orb1474920
Category	Antibodies
Description	FITC Conjugated Rabbit Anti-human IgG (H+L) secondary antibody This FITC conjugated antibody is specific for human IgG and shows no cross-reactivity with rat/mouse/goat/rabbit IgG.
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	FITC
Reactivity	Human
Concentration	1mg/ml
Immunogen	Human IgG (whole molecule).
Tested applications	FC, IF
Antibody Type	Secondary Antibody
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)