

Product Datasheet

Goat Mouse IgG (H+L) Secondary Antibody, FITC Conjugated (orb1474925)

Catalog Number	orb1474925
Category	Antibodies
Description	FITC Conjugated Goat Anti-mouse IgG (H+L) secondary antibody This FITC conjugated antibody is specific for mouse IgG and shows no cross-reactivity with human/bovine/rabbit IgG.
Clonality	Polyclonal
Species/Host	Goat
Conjugation	FITC
Reactivity	Mouse
Concentration	1mg/ml
Immunogen	Mouse IgG (whole molecule).
Tested applications	FC, IF, WB
Application notes	Western blot: 0.2-0.5 µg/mL Immunofluorescence: 10-20 µg/mL Flow Cytometry: 10-20 µg/mL
Specificity	This FITC conjugated antibody is specific for mouse IgG and shows no cross-reactivity with human/bovine/rabbit IgG.
Antibody Type	Secondary Antibody
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note

For research use only



IF analysis of human intestinal cancer tissues using mouse anti-PCK Antigen Affinity purified monoclonal antibody.FITC Conjugated Goat Anti-mouse IgG (H+L) secondary antibody was used to detect the primary antibody.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)