

Product Datasheet

Human SATB1 Protein (orb1475747)

Catalog Number	orb1475747
Description	Recombinant Human SATB1 is produced by our mammalian expression system and the target gene encoding native Met1-Asp763 is expressed with a His and a FLAG tag at the N terminus.
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-5 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	N-6His-FLAG
Note	For research use only
Uniprot ID	Q01826
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)