

## Product Datasheet

# Human Ubiquitin-conjugating enzyme E2 Z Protein (orb1475778)

<b>Catalog Number</b>	orb1475778
<b>Category</b>	Proteins
<b>Description</b>	MADE TO ORDER   Recombinant Human Ubiquitin-conjugating enzyme E2 Z is produced by E.coli expression system and the target gene encoding T218-V354 is expressed with a His tag at the C-terminus.
<b>Tag</b>	C-His
<b>UniProt ID</b>	<b>Q9H832</b>
<b>Expression System</b>	E. coli
<b>Biological Origin</b>	Human
<b>Expression Region</b>	T218-V354
<b>Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Dry Ice Shipping</b>	<b>Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.</b>
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)