

Product Datasheet

Human E3 ubiquitin-protein ligase AMFR Protein (orb1475819)

Catalog Number	orb1475819
Category	Proteins
Description	MADE TO ORDER Recombinant Human E3 ubiquitin-protein ligase AMFR is produced by E.coli expression system and the target gene encoding T450-S643 is expressed with a His tag at the C-terminus.
Tag	C-His
UniProt ID	Q9UKV5
Expression System	E. coli
Biological Origin	Human
Expression Region	T450-S643
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Dry Ice Shipping	Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)