

Product Datasheet

DWARF14 (C) Antibody, Rabbit Polyclonal (orb1478844)

Catalog Number	orb1478844
Category	Antibodies
Description	Rabbit polyclonal antibody to DWARF14 (C)
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Plant
Buffer/Preservatives	0.01% NaN ₃
Purity	affinity purification
Immunogen	C-terminal region of arabidopsis thaliana DWARF14
UniProt ID	Q9SQR3
MW	30
Tested applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558