

Product Datasheet

ABI5 (N) Antibody, Rabbit Polyclonal (orb1478895)

Catalog Number	orb1478895
Category	Antibodies
Description	Rabbit polyclonal antibody to ABI5 (N)
Clonality	Polyclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Plant
Buffer/Preservatives	0.01% NaN ₃
Purity	affinity purification
Immunogen	N-terminal region of arabidopsis thaliana ABI5
UniProt ID	Q9SJN0
MW	64
Tested applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)