

Product Datasheet

Cynomolgus CX3CL1 Protein, His Tag (orb1496211)

Catalog Number	orb1496211
Category	Proteins
Description	Cynomolgus CX3CL1 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Gly 24 - Ala 352 (Accession # A0A872TN34-1).
Target	CX3CL1
Conjugation	Unconjugated
Tag	C-10×His
Reactivity	Monkey
Form/Appearance	Powder
Buffer/Preservatives	PBS, pH 7.4
Purity	90%
UniProt ID	A0A872TN34
MW	36.5 kDa
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 36.5 kDa. The protein migrates as 50-75 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Expression System	HEK293
Endotoxins	1.0 EU per µg

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Expression Region	Gly 24 - Ala 352
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)