

Product Datasheet

NRG3 Antibody Blocking Peptide (orb1502946)

Catalog Number	orb1502946
Category	Proteins
Description	This product is NRG3 Antibody Blocking Peptide. It targets Pro-neuregulin-3, membrane-bound isoform. Purification is performed using HPLC.
Target	Pro-neuregulin-3, membrane-bound isoform
Purification	HPLC
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)