

Product Datasheet

Phospho-Smad3 (Ser423/425) Antibody Blocking Peptide (orb1505094)

Catalog Number	orb1505094
Category	Proteins
Description	This product is Phospho-Smad3 (Ser423/425) Antibody Blocking Peptide. It targets Mothers against decapentaplegic homolog 3. Purification is performed using HPLC.
Target	Mothers against decapentaplegic homolog 3
Purification	HPLC
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)