

## Product Datasheet

# phospho-SMAD3 (Ser423) Antibody Blocking Peptide (orb1506687)

<b>Catalog Number</b>	orb1506687
<b>Category</b>	Proteins
<b>Description</b>	This product is phospho-SMAD3 (Ser423) Antibody Blocking Peptide. It targets Mothers against decapentaplegic homolog 3. Purification is performed using HPLC.
<b>Target</b>	Mothers against decapentaplegic homolog 3
<b>Purification</b>	HPLC
<b>Storage</b>	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)