

Product Datasheet

TNFSF11 Antibody Blocking Peptide (orb1507324)

Catalog Number	orb1507324
Category	Proteins
Description	This product is TNFSF11 Antibody Blocking Peptide. It targets Tumor necrosis factor ligand superfamily member 11. Purification is performed using HPLC.
Target	Tumor necrosis factor ligand superfamily member 11
Purification	HPLC
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)