

Product Datasheet

SAP30L Antibody Blocking Peptide (orb1508354)

Catalog Number	orb1508354
Category	Proteins
Description	This product is SAP30L Antibody Blocking Peptide. Targets Histone deacetylase complex subunit SAP30L. Purified using HPLC.
Target	Histone deacetylase complex subunit SAP30L
Purification	HPLC
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)