

## Product Datasheet

### Catenin alpha1 Immunizing Peptide (orb15108)

<b>Catalog Number</b>	orb15108
<b>Category</b>	Proteins
<b>Description</b>	Alpha E-catenin Immunizing peptide
<b>Buffer/Preservatives</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Protein Sequence</b>	KREKQDETQTKIK
<b>Application notes</b>	This peptide is used as immunogen to produce Alpha E-catenin antibody (orb19232)
<b>Storage</b>	Shipped at ambient temperature, store at -20°C
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)