

Product Datasheet

DNAI1 Antibody Blocking Peptide (orb1515437)

Catalog Number	orb1515437
Category	Proteins
Description	This product is DNAI1 Antibody Blocking Peptide. Targets Dynein intermediate chain 1, axonemal. Purified using HPLC.
Target	Dynein intermediate chain 1, axonemal
Purification	HPLC
Storage	Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)